Genos Mini

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Internal distractor for lengthening
of the metacarpalia / metatarsalia

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Internal callus distraction for lengthening of the metacarpalia / metatarsalia

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Introduction

Thanks to Ilizarov, callus distraction has become an established procedure for lengthening and axis correction of extremities. Using an external fixator that must be left in place throughout the healing phase represents a complicated and troublesome therapy for patients. To increase patient convenience and practicability and prevent pin infections at the same time, we have developed a special mini-fixator. This internal distractor allows the use of callus distraction for correcting infantile malformations (brachymetatarsia and brachymetacarpia) and treating posttraumatic conditions. It offers high stability during the distraction and knitting phases and allows bone lengthening of up to 23 mm. Thanks to its low profile and internal application, this mini-fixator provides an excellent alternative to callus distraction by means of an external fixator.

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
Follow-up treatment

Metacarpal application

Plaster immobilization for 4–6 weeks using forearm cast with finger inclusion.

Regular x-ray checks.

From the 4th week, rehabilitation can take place on the basis of a forearm plaster shell.

Alternatively, a thermoplastic splint can be used.

Metatarsal application

Plaster immobilization for 4–6 weeks using short-leg resting cast or walking cast.

Thereafter, depending on radiological follow-up results, mobilization in postoperative shoe or ready-made shoe equipped with rigid insole.

Distraction phase

The distraction process is started after a postoperative latency phase (regenerate pre-stage) of 3–5 days. We recommend using a daily distraction length of 0.5 mm. This is achieved by two full clockwise turns of the activator carried out at different times of the day. One turn is equivalent to a distraction length of 0.25 mm. During the distraction phase, proper progress is continuously monitored by medical checks paying particular attention to the tensile forces acting on the patient's soft tissues as well as to sufficient callus formation (to be monitored by X-rays). If necessary, the distraction speed must be correspondingly adjusted.

Please note that excessive soft-tissue tension can lead to elastic deformation of the distractor. This should be monitored by regular X-ray examinations of the distraction progress. If the distraction process is carried on in spite of distractor deformation, this can have adverse consequences such as overload fracture of the device, screw fracture or screw tear-out.

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
KLS Martin do Brasil Ltda.

São Paulo · Brazil
Tel. +55 11 3554 2299
brazil@klsmartin.com

KLS Martin Medical (Shanghai) International Trading Co., Ltd

Shanghai · China
Tel. +86 21 5820 6251
info@klsmartin.com

KLS Martin India Pvt Ltd.

Chennai · India
Tel. +91 44 66 442 300
india@klsmartin.com

KLS Martin Italia S.r.l.

Milan · Italy
Tel. +39 039 605 67 31
info@klsmartin.com

KLS Martin Japan K.K.

Tokyo · Japan
Tel. +81 3 3814 1431
info@klsmartin.com

KLS Martin SE Asia Sdn. Bhd.

Penang · Malaysia
Tel. +604 261 7060
malaysia@klsmartin.com

KLS Martin de México, S.A. de C.V.

Mexico City · Mexico
Tel. +52 55 7572 0944
mexico@klsmartin.com

KLS Martin Nederland B.V.

Huizen · Netherlands
Tel. +31 35 523 45 38
infonl@klsmartin.com

KLS Martin SE & Co. KG

Moscow · Russia
Tel. +7 499 792 76 19
russia@klsmartin.com

KLS Martin Taiwan Ltd.

Taipei · Taiwan
Tel. +886 2 2325 3169
taiwan@klsmartin.com

KLS Martin SE & Co. KG

Dubai · United Arab Emirates
Tel. +971 4 454 16 55
middleeast@klsmartin.com

KLS Martin UK Ltd.

Reading · United Kingdom
Tel. +44 118 467 1500
info.uk@klsmartin.com

KLS Martin LP

Jacksonville · Florida, USA
Tel. +1 904 641 77 46
usa@klsmartin.com

KLS Martin SE Asia Sdn. Bhd.

Hanoi · Vietnam
Tel. +49 7461 706-0
info@klsmartin.com

KLS Martin SE & Co. KG

A company of the KLS Martin Group

KLS Martin Platz 1 · 78532 Tuttlingen · Germany
PO Box 60 · 78501 Tuttlingen · Germany
Tel. +49 7461 706-0 · Fax +49 7461 706-193
info@klsmartin.com · www.klsmartin.com