



## **HBS** 2 Resorb Mg and **mm.CS**

Headless Cannulated Compression Screws – resorbable



Osteosyntheses are our passion!

It is our ambition to advance these together with our clinical partners. Every day we work on developing innovative products and services which meet the highest demands on quality, and which contribute to the wellbeing of the patient.

## Table of Contents

Product Features	6 - 9
Fields of Use	10 - 11
Examples for Use	12 - 13
Product Range	14 - 31
<b>HBS 2 Resorb Mg</b>	
■ Implants	14 - 15
■ Instruments	16 - 19
■ Storage	20 - 21
<b>mm.CS</b>	
■ Implants	22 - 23
■ Instruments	24 - 29
■ Storage	30 - 31



---

## **HBS 2 Resorb Mg** and **mm.CS** Metal and yet resorbable

After over 100,000 successful HBS and HBS 2 implantations worldwide, it was time to evaluate the experience gained and start a new era – the introduction of the HBS 2 Resorb Mg compression screw. This innovative screw uses a special magnesium alloy as an absorbable material and retains the proven HBS 2 geometry so that the existing HBS 2 instruments and accessories can continue to be used.


In order to promote the advantages of magnesium technology in a targeted manner, KLS Martin is working together with Medical Magnesium GmbH, a young start-up from Aachen, which has specialized in the use of magnesium in osteosynthesis since 2018.

This collaboration offers many advantages. The resorbable magnesium implants enable optimal patient care without the need for metal removal or a second operation. In addition, these magnesium implants are used for the first time with a modified surface that slows down resorption and minimizes hydrogen gas formation, resulting in faster restoration of the original bone condition.

In addition, our users benefit from extended treatment options thanks to this collaboration. The extensive portfolio enables a wide range of procedures, particularly in foot and knee surgery. At the same time, the proven ergonomics of the HBS 2 instruments and accessories are retained, ensuring seamless integration into existing procedures.

The synergy between KLS Martin and Medical Magnesium results in an advanced material technology that combines metallic stability and resorbability in a unique innovation platform to provide novel and patient-oriented solutions in fracture stabilization.

## Feature, Function and Benefit



HBS 2 Resorb Mg and mm.CS screws are available in the following dimensions in sterile packaging:

	Total length	Length of thread, proximal	Length of thread, distal
<b>HBS 2 Resorb Mg mini</b> Ø 2.5 short thread	10 - 30 mm (1-mm steps)	2.6 mm	4 - 6 mm
<b>HBS 2 Resorb Mg mini</b> Ø 2.5 long thread	20 - 40 mm (2-mm steps)	3.7 mm	9 - 13 mm
<b>HBS 2 Resorb Mg midi</b> Ø 3.0 short thread	10 - 30 mm (1-mm steps)	3.5 mm	4 - 6 mm
<b>HBS 2 Resorb Mg midi</b> Ø 3.0 long thread	20 - 40 mm (2-mm steps)	5.0 mm	9 - 13 mm
<b>mm.CS 2.8</b>	10 - 34 mm (2-mm steps)	3.8 mm	4 - 9 mm
<b>mm.CS 3.5</b>	10 - 40 mm (2-mm steps)	4.0 mm	4 - 10 mm
<b>mm.CS 5.0</b>	30 - 60 mm (2-/5-mm steps)	4.8 mm	8 - 15 mm


The screws of all variants are made of the magnesium alloy WE43, with an additional surface modification. This leads to improved mechanical stability compared to pure magnesium or other magnesium alloys. Additionally, it allows for slower and more controlled degradation, thus providing fractures or osteotomies with sufficient mechanical support until consolidation.

During the resorption process, gas may be released, which appears as a shadow around the implant when imaging. This is a normal part of the process and should not be confused with signs of the implant becoming loose. To accurately assess the position and volume of the released gas, we recommend viewing postoperative X-rays from several planes.

## HBS 2 Resorb Mg and mm.CS – Implants


### Feature and Function

### Benefits



- Compression screw, cannulated, made of resorbable magnesium.
- In the proven HBS 2 and mm.CS geometry.
- Surgical technique for HBS 2 and HBS 2 Resorb Mg is largely identical.

- Easy placement of the screw due to cannulation.
- Can be implanted flush with the bone surface, including defined and reproducible compression.
- Higher stability compared to resorbable polymers.
- Biocompatibility according to international standards.
- Clearly visible in X-ray imaging and in CT for intra- and postoperative monitoring.
- No removal of metal/second surgery necessary.
- Continued use of HBS 2 instruments.
- Intuitive application.




- Ceramized surface (surface modification).
- Resorption process (approx. 24 to 36 months, depending on the implant size and the patient-specific metabolism).

- Retarded resorption for optimal bone healing and integration in the bone tissue.
- Minimized formation of hydrogen gas.
- Controlled and bony remodeling of the screw.
- Formation of new and compact bone material.



- Magnesium alloy (magnesium, yttrium, zirconium, rare earths).

- Rare earths as an alloy component ensure increased mechanical stability when compared to pure magnesium.



- T7/T8/T15 with self-retaining function.
- Easy pick-up, insertion and removal of the screw.

## Feature, Function and Benefit



The instruments are color-coded for easy and efficient handling. Only a few instruments are required for treatment.

Given the mechanical features of the magnesium screws, it is necessary to open the cortex prior to implantation. For this purpose, the existing HBS 2 instrument set was supplemented with a countersink. The twist drill already included in the HBS 2 set is suitable for pre-drilling as required for the magnesium screws.

This provides the option of treating indications with both HBS 2 and HBS 2 Resorb Mg using just a single instrument set.


As a user, you are also given the option of using the mm.CS magnesium implants with the corresponding mm.CS instrument set.



## HBS 2 Resorb Mg and mm.CS – Instruments and Storage

### Feature and Function


### Benefits



HBS 2 Resorb Mg mm.CS


- The overall portfolio includes three storage units:
  - HBS 2 Resorb Mg mini/midi Ø 2.5/Ø 3.0
  - mm.CS 2.8 and 3.5
  - mm.CS 5.0

- Modular design for specific use.








- Reusable countersink for opening the cortex.
- The countersink is the only addition to the existing HBS 2 set of instruments.


- Reduced torque for implanting the screws.
- Titanium and magnesium screws can be implanted using largely the same instruments and the same storage unit.



- Color-coded instruments:

- HBS 2 Resorb Mg mini Ø 2.5 
- HBS 2 Resorb Mg midi Ø 3.0 
- mm.CS 2.8 
- mm.CS 3.5 
- mm.CS 5.0 

- Easier identification of the respective instruments.



- Stainless steel storage unit in reprocessing-optimized design combined with high-performance plastic.

- High stability, but low weight.
- Large openings ensure excellent rinsing capability.
- No water residues.
- Good ergonomics.

## Step by Step to Optimal Treatment

### Fields of Use


#### **HBS** 2 Resorb Mg

The HBS 2 Resorb Mg is used for the treatment of intra-articular and extra-articular fractures and pseudarthroses of small bones and bone fragments, as well as for arthrodeses on small joints.

#### **mm.CS**

The mm.CS is used for reconstructive procedures to restore bone continuity and stabilize the correction in the osteosynthesis of fractured bones and osteotomies of small bones, for the refixation of small cortical, cancellous and osteochondral fragments or avulsion fractures and for the arthrodesis of small joints.


## Examples for Use



- ① Fractures of the processus styloideus radii
- ② Fractures of the processus styloideus ulnae
- ③ Scaphoid fractures
- ④ Metacarpal fractures
- ⑤ Phalangeal fractures
- ⑥ DIP arthrodeses



- ① Fractures of the epicondylus lateralis/medialis humeri
- ② Proximal radial head fractures



- ① Fractures of the epicondylus lateralis/medialis femoralis
- ② Bony tears of the posterior cruciate ligament
- ③ Osteochondral fractures of the patella



- ① Fracture treatment of the malleolus lateralis/medialis
- ② Tarsus/metatarsus: TMT-1 fusions
- ③ Metatarsus: Jones fractures
- ④ Metatarsus: Hallux valgus corrections with, for example, Scarf osteotomies
- ⑤ Metatarsus: Weil osteotomies
- ⑥ Phalanges: Akin osteotomy
- ⑦ Phalanges: PIP fusions

## Examples for Use in Foot Surgery: **mm.CS**

### Case study 1: Hallux valgus correction

Hallux valgus correction by means of Scarf and Akin osteotomy. Use of the mm.CS 3.5 for the distal metatarsal osteotomy and the mm.CS 2.8 in the phalanx. Imaging was performed preoperatively and at follow-up after six weeks and then three months.

Reference:

Dr. Oliver Gottschalk, Schön-Klinik Munich-Harlaching

### Case study 2: refixation of the medial malleolus osteotomy

Refixation of the medial malleolus osteotomy for treatment of an osteochondral lesion on the talus with two mm.CS 5.0, 44 mm length. Parallel ligament reconstruction with suture anchors. Imaging was performed immediately postoperatively and at the 6-week follow-up, where gas accumulation around the screws is evident.

Reference:

PD Dr. Christian Plaaß, Diakovere Annastift Hanover



Pre-op



Post-op



Six weeks post-op



Six weeks post-op



Three months post-op



Six weeks post-op

## Examples for Use in Knee Surgery: **mm.CS**

### **Case study 3: avulsion fracture of the posterior cruciate ligament**

Using a posterior approach, the avulsion fragment, which was torn out by the posterior cruciate ligament, was refixed stably with two mm.CS 3.5 of length 38 mm. The images show the fracture condition directly post surgery.

Reference:


PD Dr. Clemens Kösters, Maria-Josef-Hospital Greven

### **Case study 4: treatment of the medial femoral condyle**

Fixation of an osteochondral fragment in the femoral condyle with two mm.CS 2.8. Imaging was performed pre-op and six months post-op.

Reference:

PD Dr. Christoph Kittl, University Hospital Münster



Pre-op



Post-op




Pre-op



Six months post-op


## Implants **HBS** 2 Resorb Mg



**HBS** 2 Resorb Mg mini Ø 2.5 mm

### HBS 2 Resorb Mg mini Ø 2.5 short thread


Length of thread, proximal  
**P = 2.6 mm**




Total length (mm)	<b>STERILE</b>	Length of thread, distal D (mm)
10	53-800-10-04	4.0
11	53-800-11-04	4.0
12	53-800-12-04	4.0
13	53-800-13-04	4.0
14	53-800-14-04	5.0
15	53-800-15-04	5.0
16	53-800-16-04	5.0
17	53-800-17-04	5.0
18	53-800-18-04	5.0
19	53-800-19-04	5.0
20	53-800-20-04	5.0
21	53-800-21-04	6.0
22	53-800-22-04	6.0
23	53-800-23-04	6.0
24	53-800-24-04	6.0
25	53-800-25-04	6.0
26	53-800-26-04	6.0
27	53-800-27-04	6.0
28	53-800-28-04	6.0
29	53-800-29-04	6.0
30	53-800-30-04	6.0

### HBS 2 Resorb Mg mini Ø 2.5 long thread

Length of thread, proximal  
**P = 3.75 mm**



Total length (mm)	<b>STERILE</b>	Length of thread, distal D (mm)
20	53-801-20-04	9.0
22	53-801-22-04	9.0
24	53-801-24-04	9.0
26	53-801-26-04	10.0
28	53-801-28-04	10.0
30	53-801-30-04	11.0
32	53-801-32-04	11.0
34	53-801-34-04	12.0
36	53-801-36-04	12.0
38	53-801-38-04	13.0
40	53-801-40-04	13.0



**HBS 2 Resorb Mg midi Ø 3.0 mm**

Explanation of icons

- Magnesium
- T-Drive, cannulated
- Packaging unit

**STERILE** ↓ sterile packaged implants

**HBS 2 Resorb Mg midi Ø 3.0 short thread**

Length of thread, proximal  
**P = 3.5 mm**

Total length (mm)	<b>STERILE</b> ↓	Length of thread, distal D (mm)
10	53-802-10-04	4.0
11	53-802-11-04	4.0
12	53-802-12-04	4.0
13	53-802-13-04	4.0
14	53-802-14-04	5.0
15	53-802-15-04	5.0
16	53-802-16-04	5.0
17	53-802-17-04	5.0
18	53-802-18-04	5.0
19	53-802-19-04	5.0
20	53-802-20-04	5.0
21	53-802-21-04	6.0
22	53-802-22-04	6.0
23	53-802-23-04	6.0
24	53-802-24-04	6.0
25	53-802-25-04	6.0
26	53-802-26-04	6.0
27	53-802-27-04	6.0
28	53-802-28-04	6.0
29	53-802-29-04	6.0
30	53-802-30-04	6.0















**HBS 2 Resorb Mg midi Ø 3.0 long thread**

Length of thread, proximal  
**P = 5.0 mm**

Total length (mm)	<b>STERILE</b> ↓	Length of thread, distal D (mm)
20	53-803-20-04	9.0
22	53-803-22-04	9.0
24	53-803-24-04	9.0
26	53-803-26-04	10.0
28	53-803-28-04	10.0
30	53-803-30-04	11.0
32	53-803-32-04	11.0
34	53-803-34-04	12.0
36	53-803-36-04	12.0
38	53-803-38-04	13.0
40	53-803-40-04	13.0

## Instruments **HBS** 2 Resorb Mg

### HBS 2 Resorb Mg mini Ø 2.5 Standard Instruments:

						
1/2	1/2	1/2	1/2	1/2	1/2	1/2
26-875-02-07	26-875-13-07	26-875-00-05	26-875-06-07	26-875-17-07	26-875-03-07	53-876-00-09
K-wire guide 15 cm/6"	K-wire dispenser 15 cm/6" Ø 0.9 mm	Guide wire 125 mm Ø 0.9 mm	Depth gauge	Screwdriver T7, cannulated	Cleaning wire	Countersink, cannulated, AO attachment for HBS 2 Resorb Mg mini Ø 2.5
						





Explanation of icons

- HBS 2 Resorb Mg mini Ø 2.5
- Steel
- Silicone
- T-Drive, cannulated
- Packaging unit

**STERILE** sterile packaged implants

Twist drill, cannulated, Ø 1.9 / 0.9 mm

- 
- 
- 



1/2

Segmented step drill for DIP arthrodeses, cannulated, Ø 1.9 / 0.9 mm, AO attachment

- 
- 
- 



1/2

Guide wire double-pointed for DIP arthrodeses, Ø 0.9 mm

- 
- 
- 



1/2



1/2

Attachment	non-sterile	<b>STERILE</b>
<b>AO</b>		
with stop	26-875-19-09	
without stop	26-875-20-09	26-875-20-71
<b>cylindrical</b>		
with stop	26-875-09-09	
without stop	26-875-10-09	26-875-10-71

Length	Screw	<b>STERILE</b>
26 mm	22 mm	26-875-29-71
28 mm	24 mm	26-875-30-71
30 mm	26 mm	26-875-31-71

Length	Art. No.
60 mm	22-636-09-05
120 mm	22-637-09-05
140 mm	22-638-09-05
160 mm	22-639-09-05

26-875-01-07  
Soft tissue sleeve

- 
- 
- 
-

## Instruments **HBS** 2 Resorb Mg

### HBS 2 Resorb Mg midi Ø 3.0 Standard Instruments:



1/2

26-850-02-07

K-wire guide  
15 cm/6"



1/2

26-850-13-07

K-wire dispenser  
15 cm/6"  
Ø 1.1 mm



1/2

26-850-00-05

Guide wire  
125 mm  
Ø 1.1 mm



1/2

26-850-06-07

Depth gauge



1/2

26-850-17-07

Screwdriver T8,  
cannulated



1/2

26-850-03-07

Cleaning wire



1/2






53-876-50-09

Countersink, cannulated,  
AO attachment for  
HBS 2 Resorb Mg midi Ø 3.0





Explanation of icons

-  HBS 2 Resorb Mg midi Ø 3.0
-  Steel
-  Silicone
-  T-Drive, cannulated
-  Packaging unit

**STERILE** sterile packaged implants

Twist drill, cannulated,  
Ø 2.3 / 1.1 mm

- 
-  St
-  1 unit(s)



1/2

Attachment	non-sterile	<b>STERILE</b>
<b>AO</b>		
with stop	26-850-19-09	
without stop	26-850-20-09	26-850-20-71
<b>cylindrical</b>		
with stop	26-850-09-09	
without stop	26-850-10-09	26-850-10-71



1/2

26-850-01-07  
Soft tissue sleeve

- 
-  St
-  Sic
-  1 unit(s)

## Storage Unit **HBS 2**

The HBS 2 storage unit consists of various modules.

All HBS 2 instruments required for an operation are stored individually in the instrument tray. In addition, there is a free storage space which can be used individually.



#### Storage System

**HBS 2 Resorb Mg mini Ø 2.5**    **HBS 2 Resorb Mg midi Ø 3.0**


55-910-59-04	55-910-59-04	Lid
55-910-58-04	55-910-58-04	Insert (without contents)
55-910-56-04	55-910-56-04	Storage cage (without contents)

#### Standard Instruments


**HBS 2 Resorb Mg mini Ø 2.5**    **HBS 2 Resorb Mg midi Ø 3.0**

26-875-02-07	26-850-02-07	K-wire guide
26-875-13-07	26-850-13-07	K-wire dispenser
26-875-00-05	26-850-00-05	Guide wire
26-875-06-07	26-850-06-07	Depth gauge
26-875-01-07	26-850-01-07	Soft tissue sleeve
26-875-17-07	26-850-17-07	Screwdriver
26-875-03-07	26-850-03-07	Cleaning wire
26-875-19-09	26-850-19-09	Twist drill
53-876-00-09	53-876-50-09	Countersink

## Implants mm.CS



mm.CS 2.8



mm.CS 3.5

**mm.CS 2.8**

Length of thread, proximal  
**P = 3.8 mm**

Mg


Total length (mm)	<b>STERILE</b> ↓	Length of thread, distal D (mm)
10	53-928-10-04	4.0
12	53-928-12-04	4.0
14	53-928-14-04	4.0
16	53-928-16-04	4.0
18	53-928-18-04	5.0
20	53-928-20-04	5.0
22	53-928-22-04	6.0
24	53-928-24-04	6.0
26	53-928-26-04	7.0
28	53-928-28-04	7.0
30	53-928-30-04	8.0
32	53-928-32-04	8.0
34	53-928-34-04	9.0

**mm.CS 3.5**

Length of thread, proximal  
**P = 4.0 mm**

Mg

Total length (mm)	<b>STERILE</b> ↓	Length of thread, distal D (mm)
10	53-935-10-04	4.0
12	53-935-12-04	4.0
14	53-935-14-04	4.0
16	53-935-16-04	4.0
18	53-935-18-04	5.0
20	53-935-20-04	5.0
22	53-935-22-04	6.0
24	53-935-24-04	6.0
26	53-935-26-04	7.0
28	53-935-28-04	7.0
30	53-935-30-04	8.0
32	53-935-32-04	8.0
34	53-935-34-04	9.0
36	53-935-36-04	9.0
38	53-935-38-04	10.0
40	53-935-40-04	10.0



mm.CS 5.0

Explanation of icons

- Magnesium
- T-Drive, cannulated
- Packaging unit

**STERILE** ↓ sterile packaged implants

**mm.CS 5.0**

Length of thread, proximal  
**P = 5.0 mm**

Total length (mm)	<b>STERILE</b> ↓	Length of thread, distal D (mm)
30	53-950-30-04	8.0
32	53-950-32-04	8.0
34	53-950-34-04	9.0
36	53-950-36-04	9.0
38	53-950-38-04	10.0
40	53-950-40-04	10.0
42	53-950-42-04	11.0
44	53-950-44-04	11.0
46	53-950-46-04	12.0
48	53-950-48-04	12.0
50	53-950-50-04	14.0
55	53-950-55-04	14.0
60	53-950-60-04	15.0

## Instruments mm.CS

### mm.CS 2.8

Standard Instruments:



1/2

53-900-02-07

Universal sleeve for  
mm.CS 2.8/3.5



1/2

53-900-00-07

Drill sleeve for  
mm.CS 2.8



1/2

53-900-04-07

K-wire guide for  
mm.CS 2.8



1/2

53-900-03-07

Trocar for  
mm.CS 2.8/3.5



1/2

53-900-40-07

K-wire dispenser  
18 cm  
Ø 1.0 mm



1/2

53-900-06-07

K-wire  
150 mm  
Ø 1.0 mm







Explanation of icons

- mm.CS 2.8
- mm.CS 3.5
- Steel
- Silicone
- T-Drive, cannulated
- Packaging unit

**STERILE** sterile packaged implants



1/2

53-900-10-07  
Depth gauge for K-wires  
Ø 1.0/1.1 mm x 150 mm

- St
- 



1/2

53-900-11-07  
Drill, cannulated,  
Ø 2.2 mm,  
AO attachment for  
mm.CS 2.8

- St



1/2

53-900-08-07  
Countersink,  
cannulated,  
AO attachment for  
mm.CS 2.8/3.5

- St
- 



1/2

53-900-09-07  
Screwdriver blade T8,  
cannulated,  
AO attachment for  
mm.CS 2.8/3.5

- St
- 
- 



1/2

53-900-14-07  
MagDriver  
Handle QC M

- St
- 



1/2

53-900-15-07  
Cleaning wire  
Ø 1.0 - 1.1 mm

- St
-

## Instruments mm.CS

### mm.CS 3.5 Standard Instruments:



1/2

53-900-02-07  
Universal sleeve for  
mm.CS 2.8/3.5



1/2

53-900-01-07  
Drill sleeve for  
mm.CS 3.5



1/2

53-900-05-07  
K-wire guide for  
mm.CS 3.5



1/2

53-900-03-07  
Trocar for mm.CS  
2.8/3.5



1/2

53-900-41-07  
K-wire dispenser  
18 cm  
Ø 1.1 mm



1/2

53-900-07-07  
K-wire  
150 mm  
Ø 1.1 mm





Explanation of icons

- mm.CS 2.8
- mm.CS 3.5
- Steel
- Silicone
- T-Drive, cannulated
- Packaging unit

**STERILE** sterile packaged implants



53-900-10-07  
Depth gauge  
for K-wires  
Ø 1.0/1.1 mm x 150 mm

- 
- St
- 1 unit(s)
- 

53-900-12-07  
Drill, cannulated,  
Ø 2.7 mm,  
AO attachment for  
mm.CS 3.5

- 
- St
- 1 unit(s)

53-900-08-07  
Countersink,  
cannulated,  
AO attachment for  
mm.CS 2.8/3.5

- 
- St
- 1 unit(s)
- 

53-900-09-07  
Screwdriver blade T8,  
cannulated,  
AO attachment for  
mm.CS 2.8/3.5

- 
- St
- 1 unit(s)
- 
- 

53-900-14-07  
MagDriver  
Handle QC M

- 
- St
- Sic
- 1 unit(s)
- 

53-900-15-07  
Cleaning wire  
Ø 1.0 - 1.1 mm

- 
- St
- 1 unit(s)
-

## Instruments mm.CS

### mm.CS 5.0 Standard Instruments:



1/2

53-900-18-07  
Universal sleeve for  
mm.CS 5.0



1/2

53-900-17-07  
Drill sleeve for  
mm.CS 5.0



1/2

53-900-20-07  
K-wire guide for  
mm.CS 5.0



1/2

53-900-19-07  
Trocar for mm.CS 5.0



1/2

26-130-18-07  
K-wire dispenser  
18 cm  
Ø 1.6 mm



1/2

53-900-21-07  
K-wire  
150 mm  
Ø 1.6 mm





Explanation of icons

- mm.CS 5.0
- Steel
- Silicone
- T-Drive, cannulated
- Packaging unit

**STERILE** ↓ sterile packaged implants



1/2

53-900-24-07  
Depth gauge for  
K-wires  
Ø 1.6 x 150 mm



1/2

53-900-25-07  
Drill, cannulated,  
Ø 3.9 mm,  
AO attachment for  
mm.CS 5.0



1/2

53-900-23-07  
Countersink,  
cannulated,  
AO attachment. for  
mm.CS 5.0



1/2

53-900-22-07  
Screwdriver blade  
T15, cannulated,  
AO attachment for  
mm.CS 5.0



1/2

53-900-14-07  
MagDriver  
Handle QC M



1/2

53-900-26-07  
Cleaning wire  
Ø 1.6 mm



## Storage **mm.CS**

The mm.CS storage unit consists of various modules. One storage unit is available for the mm.CS 2.8 and 3.5 screws and a second for the mm.CS 5.0 screws.

All mm.CS instruments required for a surgical procedure are stored individually in the instrument tray.



Storage System			
mm.CS 2.8	mm.CS 3.5	mm.CS 5.0	
53-900-46-07	53-900-46-07	53-900-47-07	Lid
53-900-16-04	53-900-16-04	53-900-27-04	Storage cage (without contents)

Standard Instruments			
mm.CS 2.8	mm.CS 3.5	mm.CS 5.0	
53-900-02-07	53-900-02-07	53-900-18-07	Universal sleeve
53-900-00-07	53-900-01-07	53-900-17-07	Drill sleeve
53-900-04-07	53-900-05-07	53-900-20-07	K-wire guide
53-900-03-07	53-900-03-07	53-900-19-07	Trocar
53-900-40-07	53-900-41-07	26-130-18-07	K-wire dispenser
53-900-06-07	53-900-07-07	53-900-21-07	K-wire
53-900-10-07	53-900-10-07	53-900-24-07	Depth gauge for K-wires
53-900-11-07	53-900-12-07	53-900-25-07	Drill
53-900-08-07	53-900-08-07	53-900-23-07	Countersink
53-900-09-07	53-900-09-07	53-900-22-07	Screwdriver blade
53-900-14-07	53-900-14-07	53-900-14-07	MagDriver Handle QC M
53-900-15-07	53-900-15-07	53-900-26-07	Cleaning wire

## KLS Martin Group

### **KLS Martin Australia Pty Ltd.**

Sydney · Australia  
Tel. +61 2 9439 5316  
australia@klsmartin.com

### **KLS Martin do Brasil Ltda.**

São Paulo · Brazil  
Tel. +55 11 3554 2299  
brazil@klsmartin.com

### **KLS Martin Medical (Shanghai) International Trading Co., Ltd**

Shanghai · China  
Tel. +86 21 5820 6251  
info@klsmartin.com

### **KLS Martin India Pvt Ltd.**

Chennai · India  
Tel. +91 44 66 442 300  
india@klsmartin.com

### **KLS Martin Italia S.r.l.**

Milan · Italy  
Tel. +39 039 605 67 31  
info@klsmartin.com

### **KLS Martin Japan K.K.**

Tokyo · Japan  
Tel. +81 3 3814 1431  
info@klsmartin.com

### **KLS Martin SE Asia Sdn. Bhd.**

Penang · Malaysia  
Tel. +604 261 7060  
malaysia@klsmartin.com

### **KLS Martin de México, S.A. de C.V.**

Mexico City · Mexico  
Tel. +52 55 7572 0944  
mexico@klsmartin.com

### **KLS Martin Nederland B.V.**

Huizen · Netherlands  
Tel. +31 35 523 45 38  
infol@klsmartin.com

### **KLS Martin SE & Co. KG**

Moscow · Russia  
Tel. +7 499 792 76 19  
russia@klsmartin.com

### **KLS Martin Taiwan Ltd.**

Taipei · Taiwan  
Tel. +886 2 2325 3169  
taiwan@klsmartin.com

### **KLS Martin SE & Co. KG**

Dubai · United Arab Emirates  
Tel. +971 4 454 16 55  
middleeast@klsmartin.com

### **KLS Martin UK Ltd.**

Reading · United Kingdom  
Tel. +44 118 467 1500  
info.uk@klsmartin.com

### **KLS Martin LP**

Jacksonville · Florida, USA  
Tel. +1 904 641 77 46  
usa@klsmartin.com

### **KLS Martin SE Asia Sdn. Bhd.**

Hanoi · Vietnam  
Tel. +49 7461 706-0  
info@klsmartin.com

### **KLS Martin SE & Co. KG**

#### **A company of the KLS Martin Group**

KLS Martin Platz 1 · 78532 Tuttlingen · Germany  
PO Box 60 · 78501 Tuttlingen · Germany  
Tel. +49 7461 706-0 · Fax +49 7461 706-193  
info@klsmartin.com · www.klsmartin.com

### **Medical Magnesium GmbH**

Philipsstraße 8, 52068 Aachen, Germany  
Tel. +49 241 927845023  
info@medical-magnesium.com  
www.medical-magnesium.com